

Omnicell Pharmacy Catalogue

2021/22

The UK's
No.1 choice
for medication
adherence
packaging



Innovative, automated medication and supplies management solutions

Here at Omnicell, we provide innovative, automated medication and supplies management solutions to NHS hospitals, private hospitals, care homes and pharmacies both in the UK and globally. Our solutions help to reduce medication errors, improve patient safety and operational efficiencies so that clinical and care staff can spend more time on face to face patient care.

We are a total solution provider and our solutions range from custom-built robotic systems and automated dispensing cabinets to medication adherence packaging, automated filling machines and an eMAR system.

We offer a 'total' solution available across a variety of healthcare settings, from custom-built robotic systems and automated dispensing cabinets to medication adherence packaging, filling machines and Omnicell eMAR Solutions.

This catalogue has been specifically developed to showcase our medication adherence packaging, solutions and technology range.

Our products help improve adherence for hospitals, pharmacies, care homes, domiciliary care workers and patients.

Hospital



Central Pharmacy



Community



Retail Pharmacy



Care Home



Our portfolio crosses the entire continuum of care, from hospital to home.

Medication Adherence

A full suite of scalable medication adherence automation and packaging solutions to help pharmacies and care homes improve patient outcomes and achieve business goals.

Central Pharmacy Dispensing

Advanced robotics and intelligence to support error-free medication management and free pharmacy time for higher-value clinical work.

Intelligence

Predictive analytics, powerful insights and advisory services to help you optimize inventory, manage diversion and elevate pharmacy performance.

Point of Care

Workflow automation tools to ensure clinicians have safe, efficient and secure access to the right medication and medical supplies at the right time across all points of care.

How to order:

+44 (0)161 413 5333
ordersuk@omnicell.com
www.omnicell.co.uk

Please note, not all of our products feature in our catalogue. If you can't find what you're looking for simply contact us and our friendly customer service team will be more than happy to help.

DON'T FORGET TO FOLLOW
OMNICELL UK ON SOCIAL MEDIA



Contents

SureMed packaging

Pill packs and care packs – our card range	9
Solutions for the card range: Blisters	12
Accessories for our card range	16
Custom cards	18
Pill packs and care packs – our clear range	20
Monthly care pack kit	22
Pill pack and care pack solutions	24
Stickers and labels	28
Assembling your pack	30

Omniceil technology

EasySeal	35
DB500	36
DB Compact	37
Omniceil VBM	38
Omniceil eMAR Solution	42

Robotic & automated dispensing systems

Medimat	48
Automated dispensing cabinets	50
Terms and conditions	52
We're here to help	55



SureMed by Omnicell – quality you can trust.

- Our card and clear products are rigorously tested and carry the Class B quality standard.*
- This Class B standard means that Omnicell provides the safest packaging when it comes to the integrity and protection of repackaged drugs. It also guarantees SureMed by Omnicell are the lowest moisture permeable packs on the UK market.
- Without the right packaging, drugs can be damaged in transit or exposed to contaminants like moisture that limit their effectiveness.
- Don't take the risk – using cheaper medication adherence (MA) packaging could have worrying consequences for your patients' health and wellbeing.

SureMed by Omnicell

card & clear packaging

- Weekly pill packs
- Monthly care packs
- Blisters

*Excludes SureMed by Omnicell weekly pill pack clear basic.



SureMed by Omnicell adherence packs will help your patients take the right pills and right dose at the right time. It will also save you time and money. Our SureMed by Omnicell packs include a range of sizes, so you can pack larger quantities of medication.

Our weekly pill packs are available in both a clear and card range so you can choose the solution that is best for you and your patients. We pride ourselves on the fact that our solutions make taking medicine easy. Each personal pack contains the right pills and dose, ready to take at the right time. Contact us today and we'll support you in your discussions with your customers.



Inspired by Care: We see ways to save the planet everyday

Omnicell is a socially and environmentally friendly company.

Weekly pill pack card range & plastic template

The card in this product is sourced from legal and responsible sources.

The plastic blister is made from recycled materials.

The plastic template is made from recyclable material and can be fully recycled.

The recycled plastic used in our products is fully food grade with BRC accreditation for all suppliers of our plastic products.

Weekly pill pack clear range

The plastic tray is made from recycled materials.

The plastic tray and paper insert is full recyclable

To find out more about visit www.omnicell.co.uk

Pill packs and care packs

Our card range

We provide a variety of different pill packs in our card range for pharmacies to offer. Our weekly pill packs are designed specifically for community patients and our monthly care packs are for use in care homes. We also offer our popular weekly flexi pack with simple tear and go doses, our weekly max pack for patients with more complex medication regimes and our monthly care pack book which is split into seven days of single dose over a four week period – ideal for care home patients. Whatever your requirements, we have a design to suit you and your patients. You can even put your own brand on our packs to promote your pharmacy and services to your customers. All our packs are excellent quality, easy to use and disposable.



weekly pill pack

Our most popular 7 day, 4 dose pill pack.
pressure-seal

Qty: 250 units

Code: 3301-26A PACK



weekly pill pack max

pressure-seal

A large version of the weekly pill pack for patients taking a large number of medications

Qty: 100 units

Code: 3301-45A PACK

Please note: Our entire range of card based weekly pill packs (excluding weekly pill pack flip) are supplied complete with card, blister and filling template unless otherwise stated.

Weekly pill packs

Card range

weekly pill pack *flexi*

pressure-seal (left to right)

Allows patients or carers to tear off individual doses and carry either a single dose, or doses for the whole day.

Qty: 250 units

Code: 3315-10A PACK



weekly pill pack *flexi*

pressure-seal (right to left)

Qty: 250 units

Code: 3315-11A PACK




weekly pill pack *duo*

pressure-seal

A pack specifically designed for patients taking medication twice a day.

Qty: 250 units

Code: 3301-50A PACK



weekly pill pack *flip*

Flip style assembly.

pressure-seal (requires blisters and templates)

Qty: 250 units

Code: 3300-06 PACK

Monthly care packs

Card range

Our card range of 28 day, single-dose packs designed specifically for use in care homes.

monthly care pack *book*

(requires blisters and templates)

Book style assembly.

heat seal

Qty: 500 units

Code: 3301-01

pressure seal


Qty: 250 units

Code: 3300-02

pressure seal (bottom label)

Qty: 250 units

Code: 3300-12



monthly care pack *flip*

(requires blisters and templates)

Flip style assembly.

pressure seal

Qty: 250 units

Code: 3300-07



Solutions for the card range

Blisters

The majority of our pill packs come with blisters already included but they need to be ordered separately for use with the **weekly pill pack flip**, **monthly care pack book** and **monthly care pack flip**. These are available in a range of sizes.

weekly pill pack *flip* blisters
– small

Qty: 500 units Code: 200-47-MP




weekly pill pack *flip* blisters
– standard

Qty: 500 units Code: 200-48-MP




weekly pill pack *flip* blisters
– large

Qty: 500 units Code: 200-49-MP




weekly pill pack *flip* blisters
– extra large

Qty: 500 units Code: 200-49J-MP




monthly care pack *book* blisters
– small

Qty: 1000 units Code: 200-15-MP




monthly care pack *book* blisters
– standard

Qty: 1000 units Code: 200-16-MP




monthly care pack *book* blisters
– large

Qty: 1000 units Code: 200-18-MP



monthly care pack *book* blisters
– extra large

Qty: 1000 units Code: 200-50-MP



monthly care pack *flip* blisters
– small

Qty: 500 units Code: 200-47-MP




monthly care pack *flip* blisters
– standard

Qty: 500 units Code: 200-48-MP




monthly care pack *flip* blisters
– large

Qty: 500 units Code: 200-49-MP



monthly care pack *flip* blisters
– extra large

Qty: 500 units Code: 200-49J-MP




Extra Blisters

The majority of our pill packs come with blisters already included. You can purchase extra blisters for all of our card packs.

extra weekly pill pack blisters – standard

Qty: 250 units


Code: 200-46-MP



extra weekly pill pack blisters – large

Qty: 250 units


Code: 200-56-MP



extra weekly pill pack
max blisters

Qty: 250 units


Code: 200-65-MP



extra weekly pill pack
duo blisters

Qty: 250 units


Code: 200-53-MP



extra weekly pill pack
flexi blisters

Qty: 250 units

Code: 200-28-MP



Key features and benefits

- ✓ Available in a range of different sizes
- ✓ Environmentally friendly material
- ✓ Comply with EU legal regulations
- ✓ Produced by a British retail consortium accredited supplier

If you need advice on which blisters are right for you call us on +44 (0)161 413 5333.

Accessories for our card range

The following products should be used with our **card range** of weekly pill packs and monthly care packs.

To make sure you get the best from our pressure seal weekly pill packs and monthly care packs you will need to use a template and roller. Each box of weekly pill packs has a hard filling template inside which is fully recyclable. Rollers can be ordered and you can also order additional templates. To use the heat-seal pill pack you will need a paddle and an easy seal machine (see our technology section).



Templates

flexi

Hard Platen template for the pressure-seal weekly pill pack flexi.

Qty: 1 unit

Code: 584-69-MP

weekly pill pack

Hard Platen template for our pressure-seal weekly pill packs (excluding flexi).

Qty: 1 unit

Code: 584-64-MP

monthly care pack book

Hard Platen template for the pressure-seal monthly care pack book.

Qty: 1 unit

Code: 584-66-MP

monthly care pack book

Cork template for the pressure-seal monthly care pack book.

Qty: 1 unit

Code: 655-30

flip

Hard Platen template for the pressure-seal monthly care pack flip.

Qty: 1 unit

Code: 584-65-MP

duo

Hard Platen template for the pressure-seal weekly duo pack.

Qty: 1 unit

Code: 584-67-MP

max

Hard Platen template for the pressure-seal weekly pill pack max.

Qty: 1 unit

Code: 584-68-MP

heat-seal paddle

Qty: 1 unit

Code: 651-31

paddle book

Paddle for heat-seal monthly care pack book.



rollers

Qty: 1 unit

Code: 660-00

large roller

Roller for all card pressure-seal pill packs.



Dividers

These dividers should be used with the SureMed by Omnicell monthly care pack **card range**. They are available in both a standard and long size and a range of colours.

standard

Qty: 25 units

Code: 670-05

long

Qty: 25 units

Code: 670-06

and then add: OR for orange, YW for yellow, GR for green, W for white, RD for red, BL for blue, PK for pink

Pill pack patches

Adhesive pill pack labels. Can be used to change individual doses when patients alter medication.

Qty: 100 per roll

Code: PLUSPAK-PATCHES



Label flaps

Additional label flaps that can be attached inside the pill packs to allow additional label space.

Qty: 100 unit

Code: 720-LBLFLAPS-EB



Custom cards

We understand how important brand and reputation are, which is why we can bespoke our quality packs and print* your brand on the front cover to promote your pharmacy. We'll talk to you about what you need, create the artwork and then print them exactly as you want them.



Safety and reliability

Our card packs are second to none when it comes to the quality of the seal. We glue both sides of the card which means the seal is extremely secure and consistent.

*minimum order 25,000.

Weekly pill packs & monthly care packs

Our clear range

Our clear range of pill packs are made up of a transparent blister tray and come with seals included.

weekly pill pack *clear*

With a clear lid and tray

pressure-seal

Qty: 250 units

Code: CL01

spare seals

Qty: 20 units

Code: CL04

spare backing sheets

Qty: 100 units

Code: CL05

clear coloured seals

Qty: 250 units

Code: CL01CS

clear spare coloured seals

Qty: 20 units

Code: CL04CS

weekly pill pack *clear lite*

A transparent tray with card lid

pressure-seal

Qty: 250 units

Code: CL08

spare seals

Qty: 20 units

Code: CL09

spare backing sheets

Qty: 100 units

Code: CL10



weekly pill pack *clear max*

Large version of the weekly pill pack clear for patients with more complex medication regimes.

Qty: 100 units

Code: XL01

spare seals

Qty: 100 units

Code: XL01-S

spare backing sheets

Qty: 100 units

Code: XL05



weekly pill pack *clear duo*

A compact disposable pack for patients on twice daily medication.

Qty: 100 units

Code: N55D

spare seals

Qty: 100 units

Code: N85D

spare backing sheets

Qty: 100 units

Code: N46D



weekly pill pack *clear basic*

Our low cost pill pack. A seven day, four dose pill pack with plastic lid.

Qty: 250 units

Code: CL01-B



extra clear stickers

Can be used with any of the products from the clear range.

Qty: 250 units

Code: CL03

Monthly care pack clear kit

Our best-selling monthly care pack. The first choice for UK care homes.

monthly care pack clear

Our most popular 28 day pack. A clear tray with integrated blisters including coloured seals.

blisters and seals

Qty: 500 units

Code: CN01 and then add: XS, S, M, L, XL

additional blisters

Qty: 500 units

Code: CN99 and then add: XS, S, M, L, XL

seals only

Qty: 100 units

Code: CN98 and then add: Y, O, B, P

small roller

Qty: 1 unit

Code: N27

antibiotic seal (green)

Qty: 50 units

Code: CNG02

PRN seal (red)

Qty: 50 units

Code: CNR02

spare seals (blue)

Qty: 20 units

Code: CN02



Monthly care pack



Small roller for monthly care pack

dividers

Available in a range of colours

Qty: 50 units

Code: CN04 and then add:

P for pink
R for red
Y for yellow
O for orange
G for green
B for blue
W for white



Monthly care pack dividers

hangers

large hanger

Qty: 1 unit

Code: CN07



Monthly care pack hanger

platen

Qty: 1 unit

Code: CN08

Please note: If you are ordering the monthly care pack clear then please make sure you only order these dividers and metal hangers to use with it.

Pill pack and care pack solutions

Pilbob®

Allows easy removal of medication from pill packs.
Great solutions for patients with limited dexterity.

Qty: 25 unit

Code: 660-10



Using a Pilbob®

1. Hold the card upright.
2. Move the Pilbob from corner to corner following the edge of the blister.
3. All the medication will fall into the Pilbob and can now be taken.



Reminder cards:

For medication that does not fit on a hanger, the reminder card can be placed on the hangers for ANY of our monthly care packs to ensure a patient's dose is not forgotten.

Qty: 50 unit

Code: CN06

Alert cards:

Can be placed on a hanger for ANY of our monthly care packs where a pack has not yet been provided by the pharmacy if the medication was out of stock.

Qty: 50 unit

Code: CN18

Patient profile dividers

For use in the MAR binders.

Qty: 20 units

Code: N98

MAR sheets

Three part MAR sheets to ensure accurate communication about patient and medication regimes.

Qty: 750 units

Code: L-MAR-3

Single blank cassette medication card

Qty: 750 units

Code: N13-9007C7



Binders

Large ring binder to hold MAR sheets.

Qty: 1 unit

Code: BINDERS-THK

Standard ring binder to hold MAR sheets.

Qty: 1 unit

Code: BINDERS

Patient profile sticker

Patient labels which can be affixed to the dividers on ANY of our monthly care packs.

Qty: 250 units

Code: CN05

Patient medication time stickers

Can be affixed to top divider used with monthly care packs.

Qty: 50 units

Code: CN19

Pill pack and care pack solutions

Product inventory cards

Track product inventory. 1 Pad = 1400 labels

Qty: 1

Code: 670-01



Medication update pad

Qty: 1

Code: N39-A



Paper medication cups

10Z paper souffle pots can be used when administering medication.

Qty: 5000

Code: 670-03

Tag bag

Small, clear resealable bag.

Qty: 100

Code: N40

Tag bag label

Labels for resealable tag bags.

Qty: 50

Code: N41

Tag bag returns form pad

Qty: 1

Code: N96

Tag bag returns bundle

Pack includes: Tag bags, Tag bag label and tag bag returns form pad.

Qty: 1

Code: N40-A



White tag seals

Qty: 100

Code: N66



Black security tags

Qty: 100

Code: N66-A

Stickers and labels

On leave sticker



Qty: 250

Code: CN11

Nausea sticker



Qty: 250

Code: CN13

Hospital sticker



Qty: 250

Code: CN14

Asleep sticker



Qty: 250

Code: CN15

Refused dose sticker



Qty: 250

Code: CN09

Clinical reasons sticker



Qty: 250

Code: CN12


Sticker pack

Includes the following stickers:

- On leave
- Nausea
- Hospital
- Asleep
- Refused dose
- Clinical reasons

Qty: 6x250

Code: CN16



Blank Thermal Labels

1 roll = 3000 stickers

Size 1: 6"x8"

Qty: 1

Code: L-201

Size 2: 4"x8"

Qty: 1

Code: L-026



Max pack patches




Qty: 2 rolls

Code: MAX-PATCH

Security seals

Qty: 100


Code: N31



Cassette owner sticker

Qty: 100

Code: N38



Coloured spot stickers

Qty: 1

Code: N42



Weekly labels



Qty: 250

Code: L-902


Daily labels



Qty: 250

Code: L-903

Cassette medication card



Qty: 500

Code: N12-CAS1

Unprinted dose interval sticker

Blank dose interval stickers.

Qty: 100


Code: N37



Community patient profile labels

Qty: 100

Code: N64



Tray ID label

Qty: 100

Code: N61



Assembling your pack

Assembling the clear packs



1: Insert

Select a new empty blister tray. Stick the printed SureMed insert for pharmacy labels to the blister in the space provided in the clear tray lid. Place the correct patient medication labels on the instruction flap. Indicate the shape, colour etc. of medication under each label that is applied.



2: Dispensing

Select the appropriate medication, dispense the medication into the blister(s). For one tablet place this into the blister, if it is for two tablets then place both into the blister etc. When the blister is filled with the correct dosage of medication, the seal can be fixed to the blister.



3: Seal

Carefully peel the backing from the seal, starting from the top right hand corner. Position the seal over the cavities and lower into place. Gently smooth down with a roller to ensure the correct seal.



4: Close

When the cold seal is applied to the blister, this provides a moisture resistant seal. Be careful not to pierce any of the seal cavities. Close the clear tray lid and the pack is now complete.

- Do** store the SureMed packs and blisters in their original boxes.
- Do** store SureMed packs in a dry, indoor and air-conditioned environment.
- Do** apply firm pressure to pressure seal items.
- Do** advise the patient how to handle the pack.

- Don't** allow the SureMed packs to become damp or wet before, during or after assembly.
- Don't** fold or crease the SureMed packs except where intended along the score lines.

Storage advice


We want to make sure our products perform perfectly every time. We always recommend that our SureMed packs should be stored in their original packaging and in a dry environment. Storage temperatures should not exceed the range - 5°C to 37.5°C, be stored in direct sunlight, and should be kept away from sources of heat, and ignition. If stored properly, our products are expected to provide a quality performance for 24 months from date of production. Production dates are clearly printed on each box. It's important to rotate stock and check storage conditions meet this guidance. Please allow packs to reach room temperature before assembly (approx. 21°C).

Assembling the card packs



1: Peel

Remove the pressure sensitive liner from right to left to reveal the sealing surface. Put the card on the filling template (glue side up). The rest of the liner must be left attached to the foil side until you've completed dispensing.



2: Insert

Once the card is on the filling template, immediately insert the blister sheet into the card starting with the top row. Lower the rest of the sheet into the remaining holes. Apply the roller to ensure the blister sheet is securely held by the adhesive to the card.



3: Dispensing

Using the reverse side as a guide, fill the appropriate blister with the correct medication.



4: Remove Liner

Only when you are ready to seal the pack remove the pressure sensitive liner.



5: Fold

As soon as the liner is removed, fold the sealing surface to cover the filled blister.



6: Seal

Using the roller apply firm pressure to the whole surface to seal. Place the correct patient medication labels on the instruction flap. Indicate the shape, colour etc. of medication under each label that is applied. The SureMed by Omnicell pack is now complete.

TOP TIP:

For the perfect seal always place the filling template on a flat, even surface and remember to apply pressure using a roller!

If medication changes are needed, sealed packs can be corrected by using pill pack repair patches – see page 17.

1: Carefully pierce and remove the paper from the rear of the blister that needs to be changed.

2: Amend the medication and check as appropriate.

3: Carefully reseal the blister using a pill pack repair patch. The adhesive edge around the patch should be firmly pressed over the rear edge of the blister cavity to form a new seal. This enables medication changes to be made quickly before the pack is handed to the patient and without the need for new packs to be assembled.

Pill pack patches can be purchased from Omnicell and are available in packs of 200 patches.





SureMed by Omnicell technology

- EasySeal
- DB500 de-blistering machine
- DB Compact
- Omnicell® VBM 200F
- Omnicell eMAR Solution

At Omnicell we take health seriously, investing in technology to make life better and safer for the patient. We are passionate about improving the wellbeing of patients through pharmacy based services.

Led by a team of exceptional pharmacists, we consistently develop new and innovative solutions for our industry. This isn't just about saving time and money. It's about saving and improving the lives of your patients by providing the right pill and right dose at the right time.

It is this passion for patient care that continues to drive us forward. We've spent more than 20 years developing innovative automation to improve efficiency, patient safety and minimise the risk of hospital admission through non-adherence. Our technology saves time, money and lives.

EasySeal

The EasySeal is a revolution for heat-sealing compliance packaging.

It's pre-set digital controls, pressure indicator and electromagnetic lock; allow for a secure consistent seal to be achieved quickly and with minimal effort every time. Speeding up your pharmacy's workflow and ensuring a quality service for your customers.

Qty: 1 unit

Code: 100-70CE

- ✓ Compact design, ideal for small work spaces
- ✓ Unique auto open feature, for increased safety
- ✓ Pre-set digital controls
- ✓ Large soft grip handle for easy use
- ✓ Regulated temperature and pressure
- ✓ Consistent seal every time!



DB500 de-blistering machine

De-blister up to 280 tablets per minute with the DB500 de-blistering machine!

- ✓ The fastest manual de-blisterer available
- ✓ De-blisters the majority of tablets and capsules
- ✓ No more repetitive strain injury or sore fingers
- ✓ Quick calibration, ideal for de-blistering even small quantities
- ✓ A must for pharmacies offering adherence support to care homes or individual patients

Qty: 1 unit Code: 100-62M

Spare DB500 tray

Qty: 1 unit Code: 100-64M



DB Compact

The easy, efficient and convenient way to de-blister medication.

The DB Compact is one of the most cost effective de-blistering solutions available and its unique design means there's nothing else quite like it. Save yourself time, resource and repetitive injuries with the DB Compact de-blister machine.

- ✓ Output speed of up to 20 blisters per minute
- ✓ Frees up valuable pharmacy time
- ✓ Minimal supervision or training is required
- ✓ Improves efficiency

Qty: 1 unit

Code: 100-68M

Specification:

- Size - 300mm L x 205mm W x 120mm H
- Weight – 2.5kg
- 5 year warranty



Omnicell® VBM 200F

Automation for SureMed® adherence cards

The only automated pharmacy solution that efficiently and accurately fills and checks SureMed® adherence cards. Our filling machine ensures pharmacies have the competitive advantage to easily scale their business to improve adherence and patient outcomes today.



Compatible
with most
UK pharmacy
interfaces



Ask about our
new product,
Omnicell vMAR which
is the advanced
software integration of
two existing Omnicell
products - eMAR and
VBM.

Benefits for your business and patients

- improved accuracy and reduced checking time
- ability to track and trace each medication packed
- identifies each drug packed
- adherence cards support patients to live better for longer – improving health outcomes
- small-footprint, ergonomic design that fits through a standard doorway
- ease of access for cleaning and maintenance
- extremely intuitive and easy-to-use interface

Increase your pharmacy's revenue and decrease your expenses


- free up staff time and reduce your operating costs
- process up to 40 sealed and audited cards per hour
- assemble more packs with less resource
- supports hub and spoke model to improve efficiency and reduce costs

Why choose a pill pack?

- Clear blisters make it easy for both pharmacists and patients to see the meds, easy to check for mistakes and any mistakes can be easily amended without having to re-do the full pack.
- Customisable to your brand.
- Blisters can accommodate a maximum of 15 pills and up to 14 unique medications per card, suitable for patients on complex medication regimes.
- Fully compatible with automated machines such as the VBM 200F.
- Cards can display coloured med images and coloured time passes. Instructions can be printed on cards to enhance ease-of-use.
- Different formats of packs are available to better meet the needs of the patient.



Consolidated label with full colour drug images and dosage instructions



Ability to include patient photo as well as details for improved safety.



Code: 3301-47

SureMed by Omnicell VBM card with adhesive label

Code: 3301-48

SureMed by Omnicell VBM print direct to card



Omniceil eMAR Solution uses patient specific barcoding, coupled with original packs or blister packaging, to automate the entire care home medication management and administration process, making it **safe, simple and compliant** for everyone involved.

As well as making the medication administration process more efficient, Omnicell eMAR Solution also reduces medication administration errors, making the whole process safe, simple and compliant. CQC inspections are much easier with Omnicell eMAR Solution, providing peace of mind that all patient medication has been administered safely through easy to follow processes for all care home staff. An electronic medication audit trail verifies that everyone has received the right medication and dose at the right time, every time.

Are you looking for a safe, simple and compliant eMAR solution?

Omniceil has many benefits for your pharmacy, including:

- ✓ Generating additional revenue by investing in Omnicell eMAR Solution
- ✓ Locks care homes in for up to three years
- ✓ Integrates with your PMR or acts as a stand alone system
- ✓ Creates electronic interface with care homes - allowing the pharmacy to immediately detail any changes to prescriptions.
- ✓ Complies with NHS Digital IG Toolkit Supports compliance and CQC audit and guidelines
- ✓ The traceability and transparency of the system will improve your team's work flow - no longer waiting for care homes to locate MAR charts.
- ✓ Improve patient safety and help to provide exceptional care to your care home customers.

Omniceil provides everything you need to get started for you and your care home customers

- ✓ One easy to learn and simple to use system
- ✓ All equipment, software and licences provided
- ✓ Initial training and ongoing support every step of the way
- ✓ Pharmacy PMR interface
- ✓ 7am to 9pm technical support, 365 days a year
- ✓ Sales and marketing pack
- ✓ How to guides
- ✓ Data governance guidance to ensure NHS IG toolkit compliance
- ✓ eLearning package - to help ensure that care teams have a detailed understanding about the system and how it will help with their daily tasks (optional extra).
- ✓ Customer toolkit to help you win new business

Omniceil's eMAR Solution bundle offer can also include:

- ✓ Care home staff training
- ✓ Medication trolley
- ✓ Self medication cabinet for patients





OmniceLL eMAR Solution has many benefits for your care home customers.

OmniceLL eMAR Solution - making medication administration **safe**

Personalised resident information

- On screen photographic identification
- Resident history, allergy warnings and other information
- Full medication profile and managed changes
- Resident changes, transfers and recorded absences

Barcode scanning

- Ensures drugs and doses are correct and that medication is given to the right person at the right time, every time

Simple alert and verification system

- A reduced risk of potential medication errors with pharmacy immediately able to detail any changes to a prescription, including dose, on a centralised database which can be viewed by both pharmacy and care home.
- Missed medications or incomplete rounds
- Missed signatures
- Medication notifications
- Controlled drugs
- Irregular dosing can be managed for example Warfarin or Alendronic acid

OmniceLL eMAR Solution - making medication administration **simple**

OmniceLL provides everything your care home needs to get started

- One easy to learn and simple to use system
- All equipment, software and licences provided
- Initial training and ongoing support
- 7am to 9pm technical support, 365 days a year
- How to guides
- Data governance guidance to ensure NHS IG toolkit compliance

OmniceLL eMAR Solution - making medication administration **compliant**

OmniceLL eMAR Solution gives managers a complete picture of medication administration in their home.

Reporting suite means care homes are always ready for a CQC inspection

- Real time information and reports at the touch of a button
- Real time access to medication information from your pharmacy
- Reduces risk with medication errors
- No more paperwork or files of MAR sheets
- No illegible handwriting to try and understand
- Provides an audit trail for managers and accountability for all staff

Manager's dashboard

- Check the progress of rounds as they happen
- Track medication administration by nurse, carer or resident
- Spot missed medications, signatures and other errors fast and deal with them quickly
- Manage stock and inventory of medications in real-time
- Print off up to date eMAR sheets for urgent hospital admissions or discharge
- Managers feel in control of medication administration in their care home
- Allows managers to monitor staff performance and identify ways of improving their care home service

OmniceLL eMAR Solution helps care staff spend more quality time with residents thanks to

- Time taken for medication check-in process reduced by up to 50% (OmniceLL data)
- More efficient medication rounds
- Better medication stock control
- Up to the minute information

WATCH THE VIDEO ONLINE
www.omnicell.co.uk/emar

SureMed by Omnicell

robotic & automated dispensing systems

Medimat

The Next Generation Robot

The Medimat is Omnicell's next generation, pack picking robot that can automate the storing, management and dispensing of medication packs. Omnicell's experience of providing cutting edge technology to support retail pharmacy processes has enabled us to create a new generation of robot, one that is fit for purpose, supports effective and efficient workflows and saves you time and money.

Our next generation Medimat boasts many new cutting-edge features and functions that will revolutionise and improve how your pharmacy operates

Fully FMD Compliant



Here's what you can expect

Fill-in-box

- Automatic separation and transfer to the commissioning system.
- Fully automated and fills up to 250 packs per hour

Medimat

- Chaotic automated pharmacy storage system with internal integrated scanner.
- Modular construction allows the Medimat to fit in both smaller and larger pharmacies.
- Scan station is integrated and includes 4 cameras in the standard solution - medication packs are scanned faster than any other robot in the UK market

Dispensing Technology

Multiple models for different space layouts:

- Conveyor belt for same floor storage areas
- Pneumatic tubes for large distances
- Lift to transfer between floor levels
- Spiral slide systems, downpipes and chutes

Speedbox

- High speed channel fed stocking system for fast moving items.
- Available in three sizes to suit your dispensing volumes.
- Can be used as either stand alone solution with manual restocking or fully integrated into the Medimat with automatic stocking.
- Allows for direct or indirect dispensing.
- User can either take package directly from device's output tray or sophisticated conveyor technology will ensure medication is transported to user in the sales area.

Our experts are available to discuss any requirements you may have for our bespoke automation and robotic systems. We would be happy to visit your pharmacy and explore options that could benefit you and your business, simply contact us on +44 (0) 161 413 5333 today for more information.

Visit www.omnicell.co.uk to find out more about Omnicell's next generation medimat

Automated dispensing cabinets

Our inventory and medication management systems allow you to manage and control the usage of medicines and supplies throughout hospitals.

We offer a range of advanced but simple-to-use automated dispensing systems with drug-order integration. Each inventory management and dispensing system has been designed to control what is used, reduce stock levels and save staff time whilst enabling the hospital to improve patient safety, clinical governance and patient care.



Benefits of Omnicell automated medication cabinets:

- ✓ Eliminates mistakes caused by manual steps and workarounds
- ✓ Tracks medications throughout the hospital in real time
- ✓ Prevents drug diversion
- ✓ Improves inventory control and reduces costs
- ✓ Creates more efficient workflow processes

Contact us today

+44 (0) 161 413 5333
automationsalesuk@omnicell.com
www.omnicell.co.uk

Don't hesitate to contact us if you'd like to discuss our products and services.

We'd love to hear from you.

Terms and conditions

Omnicell Ltd

STANDARD TERMS & CONDITIONS FOR SUPPLY OF GOODS

These Conditions set out the terms on which Omnicell ("we" or "us") have agreed to supply goods to the individual, firm or company to whom our quotation is addressed or whose order is accepted by us ("you"). By signing below you confirm that you have received a copy of these terms, and that you accept that they apply in the case of any order placed by you with Omnicell.

1. Definitions

"Contract" means any contract formed between you and us for the sale and supply of Goods and/or Services which shall comprise, our written specification/ quotation/order acknowledgment, these Conditions and any other terms and conditions that we may bring to your attention.

"Consumables" means any packaging, labels, bags or similar goods or materials that are intended to be used and replaced which we supply to you.

"Customised Goods" means any Consumables that are specifically made to your specification.

"Equipment" means any heat sealing, pack filling, storage, de-blistering or similar equipment which we supply to you.

"Goods" means any Consumables, Equipment or other goods and materials we supply to you

"Price" means the price specified in our quotation or, if different, in our order acknowledgement.

"Buying Group" means a membership organisation with which we have previously negotiated and agreed trading terms and prices for the supply of Goods or Services which shall apply (subject to these Conditions) for the benefit of its members; "Representative" an employee, agent, sub-contractor or authorised representative; "Services" means any installation, support, maintenance, training or other services which we supply to you.

2. Orders and Acceptance

2.1 Each order or acceptance of a quotation for Goods and/or Services by you shall be deemed to be an offer by you to buy Goods and/or Services from us subject to these Conditions.

2.2 Subject to clause 2.3, we will not be bound by any different or additional terms or conditions or any variations to any Contract unless they are expressly accepted by us in writing.

2.3 If an order or acceptance of a quotation relates in whole or in part to Equipment and/or Services the Contract may be subject to additional terms and conditions which may be supplied to you together with our

quotation.

2.4 No order or acceptance of a quotation shall be deemed to be accepted by us until we issue a written acknowledgement of the order or (if earlier) we deliver the Goods to you.

3. Delivery

3.1 We will deliver the Goods at the place agreed in writing or, if none, at your normal place of business.

3.2 We will use all reasonable endeavours to deliver: Consumables within 3 working days (Saturday, Sunday and UK bank holidays are not working days) or such other timescale as may be agreed with the customer in writing from time to time; and Equipment within the timescale indicated in our order acknowledgement. However time for delivery of the Goods shall not be of the essence and, accordingly, we accept no responsibility or liability in the event that we are unable for any reason to keep to any applicable timescale.

3.3 When delivery is to be made by instalments, each instalment shall be regarded as a separate Contract such that any defective or non-delivery in respect of any one or more instalments shall not entitle you to refuse to accept or pay for any other instalment or Contract.

3.4 We may suspend or cease delivery of Goods at our discretion if at any time you are in breach of your obligations under these Conditions or we have reason to believe that you may not be able to pay your debts as they fall due and in any such case we shall notify you in writing whereupon we shall have no further liability or obligation to you.

4. Short Deliveries & Defects

4.1 On delivery of Goods you will be asked to sign for receipt of the consignment. The quantity of any consignment of Goods as recorded by us upon despatch shall be conclusive evidence of the quantity received by you on delivery unless you can provide conclusive evidence proving the contrary.

4.2 Any discrepancies between those Goods delivered and those ordered, and any damaged Goods must be notified, to us by 5pm on the working day following the day of delivery. To do this please call our customer service line 0161 413 5333 and speak to the Sales Administrator or

leave a message with relevant details of the problem on our answer phone at weekends or out of office hours. Should both parties agree it to be appropriate, we will then send you a damage report form, which you should complete and return to us by post.

4.3 You must provide us with a reasonable opportunity to inspect any damaged consignment and investigate any short delivery.

4.4 If you fail to give notice as specified in section 4.2 then, except in respect of any defect which is not one which would be apparent on reasonable inspection, the Goods shall conclusively be presumed to comply with our written specification/ quotation/order acknowledgement and, accordingly, you shall be deemed to have accepted the delivery of the Goods in question and we shall have no liability to you with respect to that delivery (except in relation to liability for any latent defects).

4.5 If you properly reject any delivery of the Goods that do not comply with our written specification/quotation/order acknowledgement then provided that the Goods have been properly stored by you we will at our option either replace the Goods within a reasonable time or issue a credit note at the pro rata Price against any invoice raised for such Goods. We shall have no further liability in respect thereof. Provided that the defect arises solely as a result of faulty Goods delivered by us, any replacement shall be carried out free of charge.

5. Price and Payment Terms

5.1 Invoices shall be sent in relation to each Contract the day after the Goods are despatched and/or the Services are undertaken. You shall pay the Price in the manner and at the times mentioned in our order acknowledgement or, if there is no such manner or times mentioned, not later than 30 days after the end of the month in which the invoice was issued. Time for payment shall be of the essence.

5.2 The Price specified in any quotation shall be valid for a period of 30 days from the date stated on it following which we reserve the right to vary the Price and issue a revised quotation. The presentation of a revised quotation shall supersede any

previous quotation relating to those Goods and/or Services.

5.3 Unless otherwise stated all Prices are exclusive of VAT and any other applicable taxes, duties, delivery and packaging charges, all of which amounts you shall pay in addition when paying for the Goods and/or Services.

5.4 Failure by you to pay the Price or any instalment thereof in due time shall entitle us to treat such failure as a repudiation of the Contract by you and require you to make immediate payment of all monies due or becoming due to us and to recover from you damages for such breach of contract and/or at our option to charge interest in accordance with section 5.7. Furthermore, we reserve the right to withdraw any discounts that we have agreed with you and charge our usual retail price in relation to those Goods and/or Services.

5.5 You shall make all payments due in full without any deduction whether by way of set-off, counterclaim, discount or otherwise.

5.6 We reserve the right at any time to withhold credit facilities or to limit the period of credit that we grant to you.

5.7 If payment is not made in full by the due date we may charge interest on the amount outstanding at the rate of 3% above Barclays Bank plc base rate applicable from time to time until payment is made in full.

5.8 We reserve the right to demand a non-refundable deposit payment in respect of any order you place. Where such deposit is demanded, the order will not be processed and/or accepted by us until such deposit is received.

6. Buying Groups

6.1 If you are a member of a Buying Group or have previously agreed individual terms and prices with us prior to placing an order with us for Goods or Services we shall supply the Goods or the Services to you on those terms and prices.

6.2 If you are a member of more than one Buying Group we shall endeavour, without obligation or liability, to supply the Goods or the Services to you on the best terms and prices available at the time of ordering.

6.3 In the event that:

6.3.1 a discrepancy arises between any of the terms and prices of each of your Buying Groups; or

6.3.2 a discrepancy arises between any of the terms and prices of any of your Buying Groups and any of your previously agreed individual terms and prices; or

6.3.3 it is not possible at the time of placing an order to advise you or your Representative on what constitutes the best bargain for you,

we shall not accept any responsibility for the accuracy of any advice which we give to you or your Representative at the time

of placing an order with us for Goods or Services in respect of which individually agreed or Buying Group's terms and prices represent the best bargain for you.

6.4 It shall remain your or your Representative's sole responsibility at the time of placing an order for Goods or Services to bring to our attention which individually negotiated or Buying Group's terms and prices shall be applied in respect of each type of Goods or Services ordered with us.

6.5 The terms and prices which are agreed between you or your Representative and us at the time of placing an order for Goods or Services shall be binding.

6.6 Once an order for Goods or Services has been placed and unconditionally accepted by us you shall not be entitled to receive any credit from us if you later wish to change the terms and pricing to those of an alternate Buying Group or to your individually negotiated terms and pricing in respect of a particular type of Goods or Services.

7. Storage

7.1 If we are unable through circumstances beyond our control (including your refusal to accept delivery for any reason or your failure to give adequate instructions for delivery) to deliver the Goods risk in the Goods shall pass to you, the Goods shall be deemed to have been delivered and we shall be entitled to store the Goods or to arrange storage until actual delivery is made, whereupon you shall be liable for all related costs and expenses (including storage and insurance).

7.2 If you submit an order in respect of Customised Goods and we agree to delay delivery of all or part of those Customised Goods at your request:

7.2.1 risk in the entire order of Customised Goods shall pass to you on the earlier of:

a) the date on which we send notification to you that the order is complete and available for delivery; or

b) the date when the first consignment of Customised Goods in respect of the order is actually delivered to you;

7.2.2 all Customised Goods ordered shall be deemed to have been delivered on the date risk passes in accordance with clause 7.2.1, but we shall at your request store or arrange storage for all or part of Customised Goods until actual delivery is made, whereupon you shall be liable for all related costs and expenses (including storage and insurance);

7.2.3 we may defer issuing an invoice in respect of all or part of the Customised Goods ordered until actual delivery is made but we reserve the right to issue an invoice in respect of the entire order of Customised Goods at any time after the date risk passes in accordance with clause 7.2.1.

8. Warranties

8.1 You acknowledge that you rely on your own skill and judgment in determining the suitability of the Goods for any particular purpose. We give no assurances as to suitability of Goods. You are alone responsible for ensuring that the Goods are sufficient and suitable for your purpose save only in so far as a person authorised to sign on our behalf has expressly agreed in writing that you may rely on our advice.

8.2 Subject to clauses 8.3 and 8.4, the Goods and/or Services shall be in accordance with any specification agreed in writing between you and us and shall be of satisfactory quality.

8.3 We reserve the right to alter the specification of our Goods from time to time. However, we shall notify you of any such changes and shall use reasonable endeavours to offer you a suitable replacement.

8.4 The dimensions of our Equipment are only approximate and may change, we give no warranty in respect of the dimensions and accept no liability in the event that the dimensions vary from those given in our quotation. Furthermore, we are continually developing and improving our Equipment and we reserve the right to supply you Equipment of a better specification than that which you ordered.

8.5 All warranties, conditions and other terms implied by statute or common law (save for conditions implied by section 12 of the Sale of Goods Act 1979) are, to the fullest extent permitted by law, excluded.

9. Title and Risk

9.1 Unless otherwise agreed in writing between the parties, title in the Goods will pass to you from the time they arrive at the place of delivery.

9.2 Unless otherwise agreed in writing between the parties, the Goods shall be at your risk as from the time they arrive at the place of delivery.

10. Cancellation of Order & Returns

10.1 Subject to clause 10.4 to 10.6 which are applicable in respect of Customised Goods only, you may cancel an order for all other Consumables at any time prior to the agreed delivery date. However, if you cancel within 48 hours of the agreed delivery date, we reserve the right to charge you for any reasonable costs and expenses we may have incurred in respect of your order.

10.2 Subject to the agreement with Omnicell customer service department (telephone number 0161 413 5333), Consumables, but excluding Customised Goods, that have been delivered and are no longer required or which have been ordered in error may be accepted for return within seven days of the agreed delivery date. All items must be

undamaged in their original containers, and in every respect unused and intact. Items shall be returned at your expense to such address and with such documentation as may be agreed with the Omnicell customer service department. A handling charge of 10% (or £5, whichever is the greater) is payable on each item returned.

10.3 We shall not be bound by any cancellation in respect of any Equipment unless you meet any cancellation deadlines that may have been specified or, if none are specified, then unless you notify us of such cancellation within a reasonable time before the agreed delivery date (taking into account the nature of the Equipment ordered and lead times). Any cancellation, otherwise than prior to any specified deadline, shall not affect your liability for payment in respect of the Equipment concerned.

10.4 We shall not be liable if you cancel any order for Customised Goods, further to our acceptance of your order, and you shall be and shall remain liable for payment in respect of those Customised Goods.

10.5 If you place an order with us for Customised Goods and, further to our acceptance of such order, you request and we agree that we shall manufacture and maintain an agreed volume or level of Customised Goods for your anticipated future orders, with such volumes or levels being based on historic ordering patterns or forecasts for expected use of the Customised Goods, as amended from time to time, you shall at all times be and shall remain liable for payment in respect of such back up stock of the Customised Goods.

10.6 You acknowledge and agree that if the volume or level of Customised Goods that you order from us significantly decreases to the extent that the previously agreed levels of back up stock of Customised Goods are not supplied to you, through no fault of our own, within 24 months after the date of production of the Customised Goods, we shall have the right to destroy such unused back-up stock of Customised Goods and invoice you for the same, which you acknowledge will no longer conform to the warranty provided in clause 8.2 and you shall pay such invoice in accordance with these Conditions.

11. Liability

11.1 Save as set out in clause 4.5, the following provisions set out our entire financial liability (including liability for the acts or omissions of our employees, agents and sub-contractors) to you in respect of:

- any breach of these Conditions;
- any use made or resale by you of any of the Goods, or of any product made by or incorporating any of the Goods; and
- any representation, statement or tortious act or omission including negligence arising under or in connection with the Contract.

11.2 Nothing in these Conditions shall operate so as to exclude or restrict our liability for death or personal injury caused by reason of our negligence or that of our

employees or agents.

11.3 Subject to clause 11.2, we shall not be liable for:

- any defect arising from wilful damage, negligence, failure to follow care instructions, abnormal working conditions, misuse or alteration;
- any Goods that have not been paid for by the due date;
- any defect that arises from any drawing, design or specification supplied by you; or
- any pure economic loss, loss of profit, business, revenues, anticipated savings or depletion of goodwill, in each case whether direct, indirect or consequential.

11.4 Subject to clause 11.2, our total liability to you whether in contract, tort (including negligence or breach of statutory duty), breach of third party rights or otherwise howsoever arising shall not in any event exceed the Price paid by you for the Goods and/or Services delivered or, in the case of liability which relates specifically to certain Goods or Services delivered, the Price paid by you for those Goods or Services.

11.5 Termination of any Contract for whatever reason shall not affect any rights or liabilities which have already accrued at the time of termination.

12. Force Majeure

12.1 We will not in any event be responsible for any delay in or failure to perform any Contract which is caused by factors beyond our reasonable control including acts of God, acts of governmental or any other lawful authority, strike, lock-out, blockade or other industrial action, war, riot, civil commotion, flood, fire, lightning, explosion, epidemic, trade embargo or any other act or omission on the part of any third party or sub-contractor, interruption or shortage of diesel, petrol or other fuel for transport or any failure of electricity or other fuel or equipment ("Force Majeure").

12.2 In the event of Force Majeure, we will use reasonable endeavours to remove the cause and to renew performance as soon as practicable thereafter. However, in the event that the Force Majeure persists for more than one month either party may cancel the relevant Contract with no liability to the other.

13. Assignment and Entire Agreement

13.1 You may not assign any Contract in whole or in part without our prior written consent.

13.2 We may perform any Contract or part of a Contract ourselves or assign or sub-contract to any other company.

13.3 Each Contract sets out the entire agreement and understanding between us relating to the subject matter thereof and supersedes all prior discussions between us and all prior memoranda of intent or understanding and all previous documentation in relation to such subject matter. Any variation must be signed in writing by one of our properly authorised representatives.

14. General

14.1 Any waiver by either of us for any breach of, or default under, any of these Conditions shall not be effective unless confirmed by the waiving party in writing and shall not be deemed a waiver of any subsequent breach or default and shall in no way affect the other Conditions.

14.2 Each right and remedy hereunder is without prejudice to any other right or remedy which may be available pursuant to these Conditions or otherwise.

14.3 Each provision of these Conditions is distinct and severable from the others and if at any time one or more of such provisions is or becomes invalid, unlawful or unenforceable (whether wholly or to any extent), the validity, lawfulness and enforceability of the remaining provisions (or the same provision to any other extent) shall not in any way be affected or impaired.

14.4 The parties do not intend that any Contract shall be enforceable by virtue of the Contracts (Rights of Third Parties) Act 1999 by any person that is not party to it.

15. Governing Law and Disputes

Any Contract between us shall be governed by and construed in accordance with English law and we both agree to submit to the exclusive jurisdiction of the English Courts as regards any claim or matter arising in respect thereof.

16. Notice

16.1 Any notice to be given under this agreement shall be in writing and shall be deemed duly given if left at or sent by first class recorded delivery, by e-mail or by fax, if it is sent by us to your last known address, e-mail address or fax number or if it is sent by you to:

2 Omega Drive
Riverbend Technology Centre
Irlam

Manchester

M44 5GR

0161 413 5333

ordersuk@omnicell.com

For customer service department, please post or telephone to same details above or email: ordersuk@omnicell.com Opening Hours: Monday – Friday 9.00am – 5.00pm

16.2 Any notice sent by e-mail or fax shall be deemed to have been delivered on the day of its delivery or transmission or if such day is not a working day on the next following working day thereafter unless a notification of non-delivery is received. Any notice sent by recorded delivery shall be deemed to have been received seven (7) working days after despatch and in proving the fact of despatch it shall be sufficient to show that the envelope containing such notice was properly addressed stamped and posted.



We're here to help

Our dedicated customer service and sales team are committed to helping you find the right medication adherence solution that works for you, your pharmacy and your patients. Whether you want to grow your business or simply require some advice on our products, we're here to help.

At Omnicell, we feel it's important to be your committed long-term partner and once you receive your goods from us, we are there to continue to provide support for you and your pharmacy. Our friendly customer service team are available from Monday to Friday, 8am-6pm and we're available to answer any queries or orders you may have.

With a network of sales representatives across the UK, we can visit your pharmacy within hours to support your needs and explore options that could work for you. Each one of our valued customers receives a named and experienced sales representative who is there for general advice, to take your questions and to ensure you get the most out of our products.

It's easy to get in touch with us simply:

Call +44 (0) 161 413 5333

Fax +44 (0) 161 775 6658

Email ordersuk@omnicell.com

Omnicell Medimat

The Next Generation Robot

Brand new features and functionality

- Increased capacity over the full length of the machine using Omnicell's new adapted shelving system.
- New Fill In Door Direct input mode adapted in line with FMD requirements (in addition to existing mode) so you can get items in stock much faster using the new multipicking gripper.
- FlexArea shelves to make storage flexible and adaptable – move them to where you need them most.
- Additional software options are available that offer further efficiency and safety including multi tenancy, controlled drugs options and user management.
- Runs on the latest version of Windows.
- Improved gripper technology to make dispensing faster and further reliable.
- The gripper dimensions have increased so the robot can manage larger medication packs, increasing the number of lines that can be automated.
- A new, modern look and feel. Available in four different colours and 16 standard combinations to match pharmacy branding.
- Bespoke sizes and configurations to suit your needs.
- The option to upgrade with the new Xpress module - will help your pharmacy to optimise and efficiently manage both fast and slow moving product lines that require a high output

Book your free demonstration or visit Omnicell Customer Experience Centre in Manchester by calling 0161 413 5333 or emailing ordersuk@omnicell.com

