

# FlexLock

A secure, refrigerated extension for Omnicell cabinets

Designed to offer security for supplies or refrigerated medications, The FlexLock can be connected to an Omnicell cabinet or PC. The colour touch screen on your cabinet or PC will then display the status of the FlexLock - open/closed, locked/unlocked. You can attach up to 12 FlexLocks to each Omnicell cabinet and up to 199 items can be assigned to each FlexLock allowing you to scale up or down to meet the workload of your hospital.

The FlexLock is also available with a calibrated TempCheck temperature sensor - an automated storage temperature monitor which will alert staff when the storage temperature is outside the specified range ensuring a better response time. The screen can display realtime storage temperature readings through a diagnostics window.

This product promotes increased **patient safety, workflow efficiency** and **regulatory compliance**.

## Key features and benefits

- offers maximum security
- can store up to 199 medications
- combinations available to meet the needs of your Trust
- automated monitoring of storage temperature
- real time alert system for staff



0161 413 5333  
 automationsalesuk@omnicell.com  
 www.omnicell.com/uk

