

STOCK
AVAILABLE
WHEN AND
WHERE YOU
NEED IT

Introducing... Omnicell SupplyX Inventory Management System

Bringing together a range of Omnicell supply solutions to give you a live dashboard of stock across all your sites.



Omnicell SupplyX automates the inventory management process improving traceability and providing you with real-time data at the touch of the button. It's easy to use, user friendly dashboard means you can see your current supply statistics on the move and share it easily with colleagues. It will also provide you with all the patient costing data you need so that you can make decisions on how and where to spend your supplies budget, driving efficiencies across your sites.

Improved traceability and real time data at your finger tips

- Omnicell SupplyX is a web based, real time stock level information dashboard and reporting suite. It will link to your hospital reporting system to control and report on open and closed stores for both top-up and perpetual inventory management.
- This scanning and software platform automates the whole inventory management process for hospital supplies.
- Eliminates unforeseen stock outs through re-orders based on current stock levels rather than weekly or monthly order cycles allowing you to use less and buy less.
- Improves traceability of every supply moving through multiple sites and provides you with real-time data at the touch of a button.
- Ensures stock is rotated to reduce out-of-date waste including expiry tracking and automated supplier expiry alerts.
- Omnicell SupplyX enables you to track and cost everything at a patient, procedure and consultant level.

Designed to improve efficiency and ways of working across your Trust

- Patient data capture from anaesthetics to surgery to ward will help you to understand what's spent so you can identify the most economic care pathway for your patients.
- On top of driving efficiencies, it also improves patient safety in the event of a product recall thanks to individual batch tracking and frees up clinicians time. Time spent searching overburdened shelves for items can now be better spent on face to face patient care.
- Supports hub and spoke model of stocking to gain value in large volume ordering.

Compliant with the latest NHS policies

- Omnicell SupplyX has been developed specifically around GS1 and Scan4Safety requirements.
- Omnicell is a GS1 Accredited Industry Partner and certified solution provider.
- Omnicell will help you meet efficiency demands set out in the Lord Carter report.



Contact us today to find out more

0161 413 5333 | automationsalesuk@omnicell.com | www.omnicell.co.uk

Key benefits

- Monitors all reasons for item wastage.
- Supports remote stock management.
- Single line searches for batch recalls with visibility of all used or stored items.
- Minimise wastage through expiry.
- Critical stock levels and live email alerts can be configured where stock is required urgently.
- Previous orders can be quickly recovered or order templates created and saved where repeat ordering is required.
- Eliminate nursing time spent on managing and ordering stock.
- Every aspect of GTIN, GLN, GSRN and GIA/GRAI data capture and integration has been considered to provide the very best capture, reporting and outcomes.
- Compatible with all common hospital cataloguing systems.
- Omnicell has a database of GS1 compliant GTINs from the top 15 suppliers to the NHS which can be cross-referenced to your Trust catalogue to speed up implementation.
- Product GTINs can be electronically captured via the Omnicell handheld or directly into the web based patient costing module.
- Patient GS1 compliant wristbands can be scanned in a single action, minimising keying in errors.
- Positive patient identification enforced – scan the patient barcode and you'll then see their name and date of birth for comparison as a feed from the PAS/Theatre Management System.
- Omnicell can assist with assigning a unique number to each location within the Trust. These can then be stored in the Omnicell database.
- Within each store, shelf and bin location fields can also be stored for greater visibility and these can be barcode managed.

Top dashboard features

Top-up

Order stock as you go – either as part of your routine or on the fly.

Cycle count

Check your stock levels quickly and easily.

Receipt

First line receipting and receipt stock as it comes in.

Transfer

Move your stock between different locations in your hospital.

Patient costing

Track stock usage and cost per patient from theatre to cath lab to ward.

Case picking

Pick stock based on preference card or case by case basis.

Restock

Restock store rooms and credit stock levels.

Destock

Take stock out of your system with five simple steps.

Lot and serial number

Full traceability of what has been used and what is on the shelf.

Reporting

Interactive dashboard reports for real time management information.



BY ITEM
STOCK USAGE
BREAKDOWN
AVAILABLE AT
THE TOUCH OF
A BUTTON



SCAN4SAFETY
Patient, Product, Place, Process.

Contact us today to find out more

0161 413 5333 | automationsalesuk@omnicell.com | www.omnicell.co.uk

Omnicell[®]
Inspired by Care™