

# remoteRX

Reducing picking time at the cabinet so nurses can spend more time with patients

Our remoteRX application is designed to ensure maximum efficiency of our Omnicell automated cabinets. It gives nurses the ability to remotely request medication at the patient's bedside, reducing picking time at the cabinet itself and eliminating potential queues which can be caused by a number of staff trying to process their request at the same time.

remoteRX allows medication requests to be completed at a time and location that is most convenient for the nurse, increasing nursing efficiency and patient safety by reducing the risk of any interruptions that may occur during the medication withdrawal process.

The requests are automatically sent from the remoteRX application to the appropriate Omnicell cabinet/cabinets as 'remote picks'.

Once the request has been generated, the user can go to the cabinet where they will be prompted to issue the requested medications without the need to select and pick the required items again. Users will be guided to the locations in the cabinet, of all of the remotely requested items.

## Key features and benefits

- streamlines workflows
- ensures nurses don't feel rushed while picking medication
- fully integrates with all Omnicell cabinets
- further reduces the risk for error
- frees up nurses time

