

Electronic prescribing integration

Omnicell can integrate with your electronic prescribing and medicine administration system (ePMA) to further improve safety and efficiency

The Omnicell cabinet receives a live feed of the patients' prescription profiles and includes any changes or stops to prescriptions. After choosing a patient, the nurse can then see the patient's active prescriptions, inactive prescriptions, PRN's and scheduled meds. These are filtered by the time that the prescription is due.

The ePMA integration system guides nurses to only select a patient's currently prescribed medicines from the Omnicell cabinet. For ease of use and increased safety, the cabinet has guiding lights to take the user directly to the location of the medicine in the cabinet.

Omnicell can also integrate into the user interface of your ePMA system. At the administration screen, the nurse can select the drugs due for administration and simply click 'add to Omnicell queue'. At the cabinet, the nurse just logs in and selects the patient then all the queued drugs are dispensed at the touch of a button.



NURSE CAN VIEW
PATIENT'S ACTIVE
PRESCRIPTIONS

Key Benefits of Omnicell ePMA

- **improved patient safety**
Only prescribed medicines can be issued to a patient
- **fewer medicine administration errors**
More information and safety alerts on medicines to reduce risk of errors
- **fewer missed doses of medication**
Omnicell indicates when doses of medicines are due to be given
- **fewer instances of out-of-stock medicines**
Omnicell improves stock control by automatically re-ordering
- **reduced clinical staff time spent preparing and dispensing medication, allowing more time for patients**
Medicines are easy to find and available for patients
- **fewer drugs wasted**
Stock levels can be managed easily



0161 413 5333
automationsalesuk@omnicell.com
www.omnicell.co.uk

Omnicell[®]
Inspired by Care™